



GraeL PSG

# ® graeL PSG

Industry-leading amplifier optimized for Gold Standard HD PSG studies.

Over 25 years of experience in designing PSG and EEG amplifiers, including the highly successful E-Series, has culminated in GraeL PSG. Based on the flagship GraeL amplifier, the GraeL PSG has a channel subset dedicated for Sleep studies. A true-HD PSG amplifier with channels as specified for the Gold standard AASM Type 1 studies.

HD PSG

# Technical specifications



## DIMENSIONS

Grael System Component	Size L x W x H	Weight
Amplifier	240mm x 144mm x 50mm 9.45in x 5.66in x 1.97in	700 grams 24.7 ounces

## NETWORK

POE (power over Ethernet)

## NETWORK CONNECTION

Drop in cradle for quick disconnect

## ANALOGUE TO DIGITAL CONVERTER

24 bit resolution

## DATA SAVE RATE

256, 512, 1024, 2048 sps

## SAMPLE RATE

16384 sps

## PATIENT CONNECTIONS

CF as per IEC60601-1

## SAFETY

IEC60601-1

## EMC

IEC60601-1-2

## COMMUNICATION INTERFACE

**Network Type** 802.3/802.3u twisted pair  
Ethernet with auto MDIX  
Supports IPV4 with DHCP, TCP, UDP,  
ICMP and ARP protocols

## 13 REFERENTIAL CHANNELS

**Range**  $\pm 300\text{mV}$  (resolution  $< 50\text{nV}$ )  
**Frequency response** DC to 580 Hz (user adjustable)  
**Input impedance**  $> 100\text{M}$   
**CMRR**  $> 100\text{dB}$   
**Noise (0.1 to 70Hz)**  $2\mu\text{Vpp}$  typical

## 6 DIFFERENTIAL CHANNELS

**Range**  $\pm 300\text{mV}$  (resolution  $< 50\text{nV}$ )  
**Frequency response** DC to 580 Hz (user adjustable)  
**Input impedance**  $> 20\text{M}$   
**CMRR**  $> 100\text{dB}$   
**Noise (0.1 to 70Hz)**  $2\mu\text{Vpp}$  typical

## 2 RESPIRATORY INDUCTIVE PLETHYSMOGRAPHY CHANNELS

**Sensor** Inductive band  
**Frequency response** DC to 64 Hz

## PRESSURE

**Differential input** 2 male luer  
**Range**  $\pm 100\text{cmH}_2\text{O}$   
**Frequency response** DC to 580 Hz (user adjustable)  
**Accuracy**  $\pm 1\text{cmH}_2\text{O}$  at  $25^\circ\text{C}$

## POSITION

**Sensor** Tri-axial accelerometer  
**Positions** Supine, prone, left, right, upright

## MICROPHONE

**Sensor** Electret  
**Frequency response** 5 to 580 Hz (user adjustable)

## LIGHT

**Sensor** Integrated

## OXIMETER

**SPO<sub>2</sub> range** 25 to 100%  
**Heart rate range** 18 to 300 pulses per minute  
**SPO<sub>2</sub> accuracy** 70 – 100%  $\pm 2.5$  digits within  
1 standard deviation  
**Heart rate accuracy**  $\pm 1\%$   $\pm 1$  digit

## Sample rate

**(SPO<sub>2</sub>, HR, status)** 1 sample per second, up sampled  
to match the other channels

## Sample rate

**(plethysmograph)** 256 samples per second, up sampled  
to match the other channels

## ENVIRONMENTAL REQUIREMENTS

**Temperature** Ambient operating temperature range:  
 $0^\circ\text{C} - 40^\circ\text{C}$  ( $32^\circ\text{F} - 104^\circ\text{F}$ )  
Storage temperature range:  
 $-10^\circ\text{C} - 50^\circ\text{C}$  ( $-4^\circ\text{F} - 122^\circ\text{F}$ )

## Humidity

20% - 90%-relative humidity

(non-condensing)

## Altitude

Less than 3,000m

## ELECTRICAL SPECIFICATIONS

**Power Supply** All power provided by the network  
connection using Power Over Ethernet  
(either mid-span injector or POE Switch,  
per IEEE 802.3af standard)

**Operating voltage** 48 volts

**Power consumption**  $< 10$  watts

## For more information please contact:

### Compumedics Limited, Australia:

**Headquarters**  
30-40 Flockhart Street  
Abbotsford VIC 3067, Australia  
Ph: +61 3 8420 7300  
Fax: +61 3 8420 7399  
Free Call: 1800 651 751

### Compumedics USA, Limited:

5015 West WT Harris Blvd, Suite E  
Charlotte, NC 28269  
Toll Free: +1 877 717 3975  
Ph: +1 704 749 3200  
Fax: +1 704 749 3299

### Compumedics Germany GmbH:

Global HQ for Compumedics DWL  
Josef-Schüttler-Strasse 2  
D-78224 Singen, Germany  
Ph: +49 7731 79 76 9-0  
Fax: +49 7731 79 76 9-99

All specifications are subject to change without notice.  
Please contact your Compumedics representative for  
latest technical information, pricing and product  
availability. Grael, COMPUMEDICS and the Compumedics  
logo are registered trademarks of Compumedics.

[www.compumedics.com](http://www.compumedics.com)

