

**1 AASM V 2.0\* Compliant - Profusion Sleep4 meets the latest AASM scoring rules.**

**Allows expanded reporting for hypopneas including:**

- Hypopneas meeting the new AASM rule of a 3% desaturation or an arousal
- Hypopneas meeting the traditional Medicare rule of a 4% desaturation
- Hypopneas meeting your traditional lab rule
- New report codes allow reporting of AASM, Medicare, and/or traditional hypopneas for numeric or chronological bi-level pressure treatment.

**3**




**Enhanced Pageback (PSG Online)**

**Video Review in Pageback**  
 Now available so you don't miss another minute of what the patient was doing while you attended to another.

**ECGFree trace transformation.**  
 Optional add-on allows you to score during pageback without ECG artifact.

**2**



**User Interface**

**Ribbon bar**

**Totally redesigned interface.**

Software has been redesigned to use the popular intuitive ribbon bar menu user interface, replacing the traditional menu and toolbar.

**Key benefits:**


- Related functions are now organized in logical groups - making them easier to find
- No more searching through lengthy drop-down menus!

**4**

**Scoring and Reporting**

**Key Enhancements:**

- AASM V 2.0 Compliant
- MSLT nap times now user configurable
- Added analysis and reporting support for Esophageal Pressure Monitoring (Pes)
- Cyclic Alternating Pattern (CAP) scoring and reporting



*MSLT Nap Times*


**5**

**Other Features**

**Expanded Workspaces for PSG Online** allowing four workspaces per computer and the ability to return to a default workspace when needed.

**Support for high definition cameras** which are helpful when studying patients with REM behavior and seizure disorders.

**Report expanded parameters for ASV**  
**Respiratory SUM channel can now be generated off-line**



\*AASM Manual for the Scoring of Sleep and Associated Events, Version 2.0.2 (at time of printing)

To find out more about this product, contact your local Compumedics representative.

**www.compumedics.com**

**Compumedics Limited, Australia (Headquarters):**  
 30-40 Flockhart Street  
 Abbotsford VIC 3067, Australia  
 Ph: +61 3 8420 7300 Fax: +61 3 8420 7399  
 Free Call: 1800 651 751

**Compumedics USA, Inc.:**  
 6605 West WT Harris Blvd, Suite F, Charlotte, NC 28269  
 Toll Free: +1 877 717 3975  
 Ph: +1 704 749 3200 Fax: +1 704 749 3299

