


ECG artifact removal for Profusion Sleep Software



Convert this ...

To this ... with a single click of a mouse!

ECG artifact removal at the click of a mouse!

Compumedics Sleep, using proprietary signal extraction algorithm developed by the Neuroscan research division, has developed a quick and accurate method for removing ECG artifact. Simply put, **ECGfree** identifies the ECG waveform component in each affected trace and removes that component without affecting the underlying signal of interest.

Finally!

Remove the ECG artifact from your sleep recordings.

No need to compromise signal quality by linking M1 with M2.

Retain ALL the important frequency information in each channel!

The ECGfree artifact rejection may be applied to all affected traces.

The ECGfree software plug-in option is available on all Compumedics Sleep systems using Profusion Sleep.

To find out more about this product, contact your local Compumedics representative.

www.compumedics.com

Compumedics Limited, Australia:
Headquarters
 30-40 Flockhart Street
 Abbotsford VIC 3067, Australia
 Ph: +61 3 8420 7300
 Fax: +61 3 8420 7399
 Free Call: 1800 651 751

Compumedics USA, Inc.:
 5015 West WT Harris Blvd, Suite E
 Charlotte, NC 28269
 Toll Free: +1 877 717 3975
 Ph: +1 704 749 3200
 Fax: +1 704 749 3299

Compumedics Europe GmbH:
 Werdauer Strasse 1 - 3
 01069 Dresden
 Ph: +49 351 5019-7682
 Fax: +49 351 5019-7684

Compumedics France SAS:
 Rue Jean Sapidus,
 Bât Pythagore
 67400 Illkirch-Graffenstaden
 Ph: +33(0) 981 062 869
 Fax: +33(0) 970 604 963

All specifications are subject to change without notice. Please contact your Compumedics representative for latest technical information, pricing and product availability. Compumedics ECGfree and Profusion Sleep are trademarks of Compumedics Limited. COMPUMEDICS and the Compumedics logo are registered trademarks of Compumedics Limited.