



Profusion neXus 360: A Revolution in Laboratory Management.



ACCESSIBILITY



SEAMLESS
INTEGRATION



HIGH
SECURITY



WORKFLOW
EFFICIENCY

CONTACT US
FOR A FREE TRIAL -
Nexus360@compumedics.com

Introducing Compumedics Profusion neXus 360, the next generation of Profusion neXus.

Built on the proven Profusion neXus platform with more than 15 years of customer use and thousands of users, Profusion neXus 360 offers the full functionality of Profusion neXus and more, in a fully web-based interface.

ACCESSIBILITY

- > Anywhere
- > Anytime
- > Any Device
(supporting HTML5 browser)

SEAMLESS INTEGRATION

- > Seamless hardware and software with user-configurable reporting
- > Fully managed by Compumedics
- > Scalability
- > HL7 Support
- > Multi-site management

HIGH SECURITY

- > Digitally secure study "sign-off"
- > Two-factor authenticated login
- > HIPAA compatible
- > All web traffic is securely encrypted

WORKFLOW EFFICIENCY


- > Simple management of access privileges
- > Web-based review and reporting
- > Automated updates and backups
- > Dynamic scalability to suit growing labs

PROFUSION NEXUS 360 FEATURES:

- > Simple, browser/internet-based access via HTML5
- > Two-factor Authentication Access
- > Digitally secure study "sign-off"
- > User-defined, group-based access privileges
- > Template/Document Integration
- > Non-editable audit-log
- > Multi-language Support (English, French, Chinese, Spanish)
- > Fully managed Cloud Service, simple installation, reliable system backups and easy system updating
- > In-lab acquisition and real-time uploading to the web

PLATFORM AND BROWSER INDEPENDENT:

PHONE



TABLET



PC/COMPUTER



*Optimized for touch displays

FOR MORE INFORMATION PLEASE CONTACT:

Peter Parslow Ph.D
eHealth Commercial Manager - (Located in Sydney, AUSTRALIA)
Cell: +61 (0)405 053 141 Email: peter.parslow@compumedics.com

Compumedics Limited, Australia:
Headquarters
30-40 Flockhart Street
Abbotsford VIC 3067, Australia
Ph: +61 3 8420 7300
Fax: +61 3 8420 7399
Free Call: 1800 651 751

Compumedics USA, Inc.:
5015 West WT Harris Blvd, Suite E
Charlotte, NC 28269
Toll Free: +1 877 717 3975
Ph: +1 704 749 3200
Fax: +1 704 749 3299

Compumedics Europe GmbH:
Werdauer Strasse 1 - 3
01069 Dresden
Ph: +49 351 5019-7682
Fax: +49 351 5019-7684

Compumedics France SAS:
2 rue Henry Bergson,
F-67000 Strasbourg
Ph: +33(0) 981 062 869
Fax: +33(0) 972 643 875