



Full Spectrum Neurodiagnostics - from Clinical to Research





Compumedics offers a comprehensive range of solutions for the Neurology facility – from clinical epilepsy and long-term home monitoring to research EEG, Compumedics has a solution for you.

As always, Compumedics is your trusted partner in fulfilling your business and research aspirations.



Explore our Solutions

Clinical, LTM and Research Solutions



Okti®

Everything you want.
Simple, Portable, Wireless
and High-Definition.

Okti® is the world's first high-definition portable EEG amplifier with interchangeable modules that allow you to record standard clinical montage and high-density EEG - up to 128 channels with one device.



Okti Onsite® A.V.S.

EEG with video recording
in "Real Time"!

Monitor your patient's home ambulatory studies with CONFIDENCE! Access through any commercial browser via the internet, in the home, office or mobile.

Grael® Range

HD - PSG / EEG



Compumedics' Grael® Range is a technologically advanced solution for both clinical and research applications. Choose the model that best suits your needs: Grael®4K for dual platform EEG and PSG, or Grael® EEG and LT for dedicated EEG studies.



SynAmps RT®

The Gold Standard Amplifier
for high-end Research.

SynAmps RT™ has been engineered to meet the rigors of the most demanding recording environments in brain research. With 64 to 512 channels, sample rates up to 20KHz and a host of options including triggers and STIM.


MEG Solution

OrionLifeSpan™ MEG


Opening NEW FRONTIERS in Clinical and Research MEG.

A new frontier in MEG, Orion Lifespan™ is a rotating two-in-one dual helmet EEG/MEG system for both adult and pediatric applications.

Compumedics has overcome the barriers of earlier MEG technology with the introduction of the Orion LifeSpan™ MEG. It provides increased precision and accuracy, coupled with the powerful, fully-integrated CURRY® analysis software.



World Class Software



Profusion EEG™ Software

World Class EEG
Analysis Software Suite


Profusion EEG™ - Add true flexibility to your equipment with full EEG functionality - with an advanced user-friendly interface. Provides the user with reports, mapping, template matching and multiple D1 resolution video windows. Functions with all Compumedics EEG Amplifiers.



Profusion Nexus® 360

Web-based Laboratory
Management System

Profusion neXus 360™ is a web-based patient data and Lab Management System that provides integrated hardware and software solutions for both Sleep and Neurology clinics - anywhere, anytime and on any device.



CURRY® Software

Multi-modal
Neuroimaging software

The Worlds most powerful Multi-modal Neuroimaging software. CURRY® sets the standard for source localization and visualization tools, integrating structural and functional image modalities (EEG, MEG, MRI, fMRI, CT, PET & SPECT) into a single 3-D image. Easy to use and amazingly functional, it is clear why Neuroscan remains at the pinnacle of Neurological Research.



For more information contact -

Email: customerservice@compumedicsusa.com

Customer Support: support@compumedics.com

USA Toll Free: 877 717 3975



www.compumedics.com



www.compumedics.com

All specifications are subject to change without notice.
Check with your local Compumedics representative for product availability in your region.
Okti, Onsight, Neuvo, Profusion, Nexus, Graef, SynAmps, CURRY, OrionLifeSpan,
Compumedics and the Compumedics logo are all registered trademarks
of Compumedics Limited.

AL318 Issue 0 - USA